

PRODUCT INFORMATION

PVC Gown With Thumb Loop

Description:

Our lightweight soft PVC Gown with thumb loops offer versatile general purpose barrier protection making them a fantastic choice for a whole range of tasks.

- ✓ Waist back tie
- ✓ Thumb loops help secure the sleeve along arms
- ✓ Universal size
- ✓ Non sterile
- ✓ Latex free
- ✓ Overall length 122cms
- ✓ Over shoulder protection

Type:

- Soft polyvinyl chloride

Applications:

Hygiene rules, staff protection, dispensing, education, carers

Packaging: 10 pieces per bag - 10 bags per case
100 pieces case quantity

Product Codes:

41513	One Size	0.90mm - Blue
-------	----------	---------------

Size Chart:

Measurements	cm's
Width (cm)	90
Length (cm)	122

Colourways:


 Blue

Standards:



PPE Regulation (EU) 2016/425
Plastic Materials and Articles in Contact with Food Regulation (EU) No 10/2011

Outer Carton Artwork:



The information shown on this specification sheet is intended as a guide and is advisory only. All details were correct at the time of issue. As conditions of use are ultimately beyond our control, we advise that the product is tested for the particular application before use.

Issue No.: 1

Issue Date: 14.05.2020