

THE NEW & IMPROVED FIXSTEP WALKER

Our new and improved fixed walker available in tall and short sizes with high volume pump air liners or standard liners.

Open orders for the current Jura Fixed Walker will be **upgraded** to the new improved product at no extra charge!





Two-Part Injection Moulded Sole Unit
Improves shock absorption and creates less impact when walking



Strap Lock System
Prevents the straps from rotating or moving out of place when walking



Flared & Deeper Heel Section
Retains the foot and reduces pressure on Achilles area



Moulded High Volume Pump on Air Liner
Provides balanced compression on both sides of the ankle



High Quality Hook/Loop Straps
Designed for easier patient opening and closure



Reinforced Moulded Plastic Uprights
Provide increased durability and contain no metal, making them x-ray safe