

RDS¹

1mL reduced dead space syringe



Sterile 1mL reduced dead space syringe for single use.

Supplied in boxes of 100 syringes. An outer carton contains 36 boxes, totalling 3,600 syringes (see page 3 for details).

Product code

LW1

Materials

Syringe body: Medical grade Polypropylene (PP 5250T)

Plunger gasket: Isoprene rubber (latex-free)

Syringe barrel and gasket lubricating oil: Medical grade Silicone oil (DC360)

None of these materials contain recycled content

Sterilisation

Sterilised by: Ethylene oxide gas

Conforming to: ISO:11135 Sterilisation of medical devices
(validation and routine control of ethylene oxide sterilisation)

Time to expiry: 5 years from manufacture

Product conformance

Class IIa Medical Device

Fully compliant with ISO:7886-1: Sterile hypodermic syringes for single use

European Medical Device Directive 93/42/EEC (CE Class:1)

Manufacturing compliant with Quality Assurance Standards ISO:13485

Compatibility

All Luer fit needles excluding Total Dose / low dead space needles

Sterifilt compatible

Manufacture and distribution

RDS¹ 1mL reduced dead space syringes are manufactured exclusively for Exchange Supplies, UK.

1, Great Western Industrial Centre, Dorchester Dorset DT1 1RD. www.exchangesupplies.org

Manufactured for	Shanghai Kindly Enterprise Development Group Co., Ltd.
Exchange Supplies by:	No. 658 Gaochao Road, 201803 Shanghai, People's Republic of China.

EC representative:	Shanghai International Holding Corp. GmbH (Europe). Eiffestraße 80, 20537 Hamburg, Germany. T: 0049-40-2513175. F: 0049-40-255726
---------------------------	--

KDL UKCA representative:	KDL – Overseas Business Development Director. Unit 14 Cowley Mill Road, Uxbridge, Greater London UB8 2DB.
---------------------------------	--



Graduations in
0.01mL increments

RDS¹ 1mL reduced dead space syringe

3

Packaging

Item	Weight (approx.)	Material(s)
Blister	3.49g	Printed surface: Medical grade 63gsm coated paper Blister: 0.09mm thickness medical complex film
Box 100	430g	E corrugated card
Carton 3,600	14.5kg	Uncoated 630gsm AB flute double wall card. Consisting of: 150gsm 3A Kraft Test Liner printed outer layer, with 120gsm and 100gsm A grade Senior corrugated and 140gsm B grade Kraft Test inner Liner. Recyclable.

None of these materials contain recycled content.



Blister

160.5mm x 35mm
GTIN: 06930197835731



Box 100

210mm x 165mm x 90mm
GTIN: 16930197835738



Carton 3,600

580mm x 440mm x 530mm
GTIN: 26930197835735

EXCHANGE
S U P P L I E S

1 Great Western Industrial Centre, Dorchester,
Dorset DT1 1RD | T: +44 (0) 1305 262244