

PRODUCT INFORMATION

Powderfree Vinyl Gloves

Description:

Made from high grade Polyvinyl Chloride, our Powderfree Vinyl Gloves are a simple and effective glove, ideal for clinical and automotive use alike. And as they're latex free, they are perfect for latex allergy sufferers.

- Medical Grade – AQL 1.5 – conforms to EN 455 parts 1-4
- Powder free
- Latex free
- Smooth finish
- Beaded cuff

Type:

- High quality polyvinyl chloride

Applications:

Tasks with clinical tasks, laboratories, food processing, food handling, automotive, assembly work



Colourways:

C Clear


Product Codes:

11401-4 S - XL Vinyl Powder Free - Clear

Standards:



PPE Regulation (EU) 2016/425
Plastic Materials and Articles in Contact with Food Regulation (EU) No 10/2011
BS EN 455 Parts 1-4:2000 to 2009 - Medical Gloves for Single Use



Size Chart:

Measurements in mm	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)
Hand Length (A)	171	182	192	204	215
Hand Circumference (B)	178	203	229	254	279

Issue No.: 4

Issue Date: 24.04.2020

SUPERTOUCH®

Packaging: 100 pieces per dispenser - 10 dispensers per case

Outer Carton Artwork:



Feature Symbols:



Product is free from powder



Product does not contain latex



Single use only
Do not reuse

The information shown on this specification sheet is intended as a guide and is advisory only. All details were correct at the time of issue. As conditions of use are ultimately beyond our control, we advise that the product is tested for the particular application before use.