

Readigloves® Nytraguard® Bluple®

Advanced Formula Nitrile Gloves

Nytraguard® Bluple® is no ordinary glove; it provides the latest in thin wall nitrile technology without compromising on strength. An innovation in high performance nitrile gloves, this medical grade, super lightweight glove gives you high comfort, great sensitivity and toughness - an unbeatable combination and in a distinctive new colour - Bluple!

- **Super lightweight, thin wall construction** - very high sensitivity
- **High stretch** - very high degree of comfort and fit
- **High durability** - tear and puncture resistance
- **Maximum biological resistance** - for protection
- **Medium grip** - suitable for multiple tasks, fingertip textured
- **Easy donning** - for speed and convenience
- **Powder free** - for hygiene and convenience
- **Latex free** - no risk of natural rubber latex allergic reaction
- **Food contact approved** - safe for use with all food types
- **Colour** - blue/purple

Nytraguard® Bluple® is the optimum solution for high volume nitrile glove requirements in acute and community healthcare settings, but equally at home for a variety of non clinical tasks. Because of its high sensitivity and dexterity this glove is ideal for dental procedures. Nytraguard® Bluple® gloves are available in cartons of 100 as standard, or in more economical 200 glove cartons, for cost efficiency and reduced frequency of replacement.

Technical Performance	
Classification:	Single use medical examination glove, Non sterile Class I Medical Device
Materials:	Acrylonitrile butadiene; Does not contain latex Does not contain styrene; Powder free
Design:	Ambidextrous, finger textured outer surface, beaded cuff
Inspection Criteria:	AQL 1,5
Shelf Life:	3 years
Regulatory Compliance:	<p>Medical Devices Directive 93/42/EEC</p> <ul style="list-style-type: none"> - BS EN 455, Parts 1, 2, 3 and 4. Medical gloves for single use <ul style="list-style-type: none"> - Minimum force at break - 6 Newtons - EN ISO 10993, Biological evaluation of medical devices <p>Personal Protective Equipment Directive 89/686/EEC</p> <ul style="list-style-type: none"> - Simple design - BS EN 420. Protective gloves - general requirements and test methods <ul style="list-style-type: none"> - Gloves are 240mm in length for medical applications - Dexterity level 5 (high dexterity) - BS EN 374 parts 1, 2 and 3. Protective gloves against chemicals and micro organisms <ul style="list-style-type: none"> - Micro organism resistant - Low level chemical protection <p>Regulation (EC) No 1935/2004 and Regulation (EU) No 10/2011. Plastic materials and articles intended to come into contact with food</p> <ul style="list-style-type: none"> - BS EN 1186. Materials and articles in contact with foodstuffs - plastics - Suitable for all food types
Applications:	Hospital - medium to high risk clinical procedures and routine tasks, community nursing, GP surgeries, dental surgeries, emergency services, nursing homes, care homes, domiciliary care, general use, cleaning, light industrial, local authority, police, customs, education, food and catering



Product Code	Size	Gloves per Carton	Unit of Sale
9535	Extra Small	100	10 Cartons
9536	Small	100	10 Cartons
9537	Medium	100	10 Cartons
9538	Large	100	10 Cartons
9539	Extra Large	100	10 Cartons
9545	Extra Small	200	10 Cartons
9546	Small	200	10 Cartons
9547	Medium	200	10 Cartons
9548	Large	200	10 Cartons
9549	Extra Large	200	10 Cartons

