



Experts in  
Medical Refrigeration

# 353L Freestanding Pharmacy Fridge



## Key Features

- ✓ Represents your vaccine vial
- ✓ Provides you with the true temperature of your vaccines
- ✓ Ensures you always measure what matters
- ✓ Built-in warning system for safeguarding vaccines
- ✓ Now with 5 year warranty as standard \*UK only



\*UK Only



**Control**  
TRUE TEMPERATURE



## PGRC353UK Pharmacy Fridge – Glass Door

<b>SKU Code</b>	<b>444443786</b>
EAN Code	5052263037861
Gross Capacity (Litres) / Net Capacity (Litres)	353 / 349
Unit External Dimensions W / D / H (mm)	595 / 625, 660 with handle / 1860
Product Weight (kg)	79
Internal Dimensions W / D / H (mm)	480 / 445 / 1682 (step:175D x 265H)
Temperature Range	+2°C to +8°C
Secondary Probe	Yes
Bottle Contents / Size	Silicone / 50ml
External Finish	Painted Aluminium
Internal Finish	HIPS*
Type of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish / Type	Painted Aluminium / Glass
Door Style	Hinged, Right Handed only, Non-reversible
Door Open Depth (mm)	1180
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	2 Adjustable Front Feet, 2 Rear Rollers
Temperature Display on Panel / Controller Type	Yes / Digital, External
Temperature Alarm	Yes
Number of Shelves / Shelf Positions / Shelf Material	6 / 14 / Wire
Shelf Size W / D / H (mm)	385 / 500 / 6
Number of Baskets / Basket Material	1 / Wire
Interior Light / Fan / Spark free	Yes / Yes / No
Auto or Manual Defrost	Automatic
Climate Class	SN
Refrigerant Type	R600a
CFC Free / HFC Free / HCFC Free / Ammonia Free	Yes
UK Plug / Other Plug Types Available	Yes / Yes
Voltage Supply	220-240
Aux Output	Yes
Insulation Thickness (mm) / Insulation Type	60 / Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	Yes