

CE marked Face Shield for Covid-19 UK manufactured, reusable, robust and comfortable.




Used, proven and popular in hospitals, GPs and social care

Product Advantages

- Designed for comfort over many hours of use
- Designed to be cleaned and sterilised to enable reuse
- Integrated ear-savers protect against discomfort from mask straps
- Industry leading sparkling optically clear visor coated to reduce fogging and static
- Simple to adjust headband, easy to secure and fit around tied back hair
- Only 3 components, assembled at point of use or supplied ready assembled
- Manufactured in the UK from recyclable materials
- Optimised for manufacture of up to 10,000/day


Full face shield



Design features

DELIVERY

- Delivered flat to reduce package size and storage requirements
- Standard 1mx1m pallet delivers 3000 SheffShields
- Standard box sizes of 50 or 250 units
- Option to be delivered assembled as units where rapid deployment and use is needed



Certification: BSI's PPE for Healthcare Professionals 2020/403—Eye protection Technical Specification
Notified Body: BSI Assurance UK Ltd
Notified Body No.: 0086 **BSI Certificate No.:** CE729058
This protective face shield is manufactured for COVID-protection only.
This protective face shield is not a PPE device for general use and shall not be used for purposes other than protection against COVID-19.

