SUPERTOUCH















CUSTOMER INFORMATION

Powderfree Vinyl Gloves

Description:

Made from high grade Polyvinyl Chloride, our Powderfree Vinyl Gloves are a simple and effective glove, ideal for clinical and automotive use alike. And as they're latex free, they are perfect for latex allergy sufferers.

- Industrial Grade AQL 4.0
- Powder free
- Latex free
- Smooth finish
- Beaded cuff

Fabric:

- High quality polyvinyl chloride

Applications:

Tasks with clinical tasks, laboratories, food processing, food handling, automotive, assembly work

Colourways:

Clear









Green

Product Codes:


11201-4	S-XL	Vinyl Powder Free - Clear
11211-4	S-XL	Vinyl Powder Free - Blue
11231-4	S - XL	Vinyl Powder Free - Green

Standards:





EC Council Directive 89/686/EEC Plastic Materials and Articles in Contact with Food Regulation (EU) No 10/2011



Size Chart:

Measurements in mm	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)
Hand Length (A)	171	182	192	204	215
Hand Circumference (B)	178	203	229	254	279

Issue No.: 2 Issue Date: 15/05/2014






Packaging: 100 pieces per dispenser - 10 dispensers per case

Outer Carton Artwork:



Inner Artwork:



Feature Symbols:



Product is free from powder



Product does not contain latex



Single use only Do not reuse

The information shown on this specification sheet is intended as a guide and is advisory only. All details were correct at the time of issue. As conditions of use are ultimately beyond our control, we advise that the product is tested for the particular application before use.

SUPERTOUCH





CUSTOMER INFORMATION

Powderfree Vinyl Gloves

Description:

Made from high grade Polyvinyl Chloride, our Powderfree Vinyl Gloves are a simple and effective glove, ideal for clinical and automotive use alike. And as they're latex free, they are perfect for latex allergy sufferers.

- Medical Grade AQL 1.5 conforms to EN 455 parts 1-4
- Powder free
- Latex free
- Smooth finish
- Beaded cuff

Fabric:

- High quality polyvinyl chloride

Applications:

Tasks with clinical tasks, laboratories, food processing, food handling, automotive, assembly work

Colourways:

Clear

Product Codes:


11401-4 S-XL Vinyl Powder Free - Clear

Standards:





EC Council Directive 89/686/EEC Plastic Materials and Articles in Contact with Food Regulation (EU) No 10/2011 BS EN 455 Parts 1-4:2000 to 2009 - Medical Gloves for Single Use



Size Chart:

Measurements in mm	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)
Hand Length (A)	171	182	192	204	215
Hand Circumference (B)	178	203	229	254	279

Issue No.: 2 Issue Date: 15/05/2014






Packaging: 100 pieces per dispenser - 10 dispensers per case

Outer Carton Artwork:



Inner Artwork:



Feature Symbols:



Product is free from powder



Product does not contain latex



Single use only Do not reuse

The information shown on this specification sheet is intended as a guide and is advisory only. All details were correct at the time of issue. As conditions of use are ultimately beyond our control, we advise that the product is tested for the particular application before use.